

# Aromatic Compounds Analysis

## 1. Analytical Condition

### YL 6500 Series

Oven : 35°C (for 8 min)- temp-programed 10 °C/min-230 °C(for 5 min)

Column : HP-INNOWAX(30m\*0.32mm\*0.5um)

Carrier gas : He, 1.0 ml/min (Split ratio 20:1)

Injector : Capillary 250°C

Detector : FID 250°C

Sample : BTEX, Dichlorobenzene(Gas)

Injection Volume : 2ml (Gas tight syringe)

## 2. Chromatogram

